

COMPUTERISED

# AUTOMATED PERIMETER

▶ *AUTOMATED PERIMETER AP-1000 / AP-2000*

- *Modern Eye Tracking System*
- *LED-Technology*
- *Digital Camera Controlled*
- *Blue On Yellow Test*
- *Multilingual User Interface*
- *Customised Test Programs*
- *Driving Test*



# STATE of the ART

SOPHISTICATED PERIMETRY WITH THE TOMÉY AP-1000 / AP-2000



Both Toméy Perimeters are developed to meet the needs of modern Ophthalmology. The AP-1000 is designed to fulfil all standard tests and it is therefore economically extreme attractive. Using the state of the art available technology - such as LED-illumination, blue on yellow tests and electronic assisted patient positioning - the AP-2000 sets the latest standards in modern perimetry.

## *Modern Eye Tracking System*

Both Systems are equipped with a "state of the art" digital passive eye tracking system. It records only measurements taken under perfect fixation conditions. Rejected recordings will be automatically repeated. The patients eye is displayed at any time during the examination for visual control by the examiner. Of course, the fixation control with Heijl-Krakau method is also integrated.

## *Short Examination Time*

Due to the fact that the units can take former testing protocols of the same patient into consideration, the examination time can be reduced tremendously. Already known defects will be screened and not necessarily checked through a complete testing routine again.


## *Pupil Measurement*

You can measure the pupil diameter at any time even before, after or during the examination.


## *Multilingual User Interface*

Both Perimeters are equipped with a multilingual user interface. Additional languages can be easily integrated.





3D-ANALYSIS AND CROSS SECTION



TOMEY PERIMETER PRINTOUT

## Customised Test Programs

Beside the standard tests you have the possibility to define new customised examination programs or change the parameters for your individual needs. An unlimited number of customised test programs can be added.

## Driving Test

The integrated standard driving test can be modified according to your individual demands.

## Modern Data Handling


The AP-1000 comes with a modern cube PC that communicates with all USB printers. All acquired data can be transmitted to a network and sent in many different formats depending on your requirements. The design of the AP-2000 is more compact due to the integrated PC. This saves valuable space and reduces cables.

## LED-Technology / Blue On Yellow Test


The AP-2000 integrates the LED background illumination and provides, using Goldman V fixation spot size, a standardized blue on yellow test routine.

## Electronic Chin Rest

The AP-2000 is equipped with an electronic chin rest for exact patient positioning.



Single Map Screenshot



Visible Difference Map (printout)

## AP-1000 AUTOMATED PERIMETER

## AP-2000 AUTOMATED PERIMETER

### SPECIFICATIONS

<b>Stimulator Bowl Type</b>	Part hemispherical, radius 300 mm.
<b>Visual Field Extent</b>	Integrated diffusing surface 100°

#### Standard Test Fields

<b>Full 50°</b>	164 points
<b>Glaucoma 22°/50°</b>	104 points
<b>Central 30°</b>	120 points
<b>Central 22°</b>	96 points
<b>Wide 22°/30°</b>	128 points
<b>Peripheral 30° to 50°</b>	72 points
<b>Macula 10°</b>	48 points
<b>Driving 50°/80°</b>	192 points

<b>Stimulus Source</b>	Rear projection LED
<b>Stimulus Colour</b>	Green 570 nm
<b>Stimulus Size</b>	Goldmann size III (0.43°)
<b>Stimulus Intensity</b>	0.03 asb to 1000 asb in 15x3 dB or 45x1 dB steps
<b>Exposure Time</b>	Adjustable: 0.1 to 9.9 s
<b>Response Time</b>	Adjustable: 0.1 to 9.9 s
<b>Inter Test Delay</b>	Adjustable: 0.1 to 9.9 s
<b>Background Illumination</b>	10 asb (3.2 cd/m), automatic level control
<b>Fixation Control Method</b>	Heijl-Krakau - blind spot method and digital camera
<b>Fixation Monitor</b>	Yellow LED, 588 nm
<b>Test Lens Diameter</b>	38 mm

#### Dimensions & Electric Requirements

<b>Dimensions WDH</b>	740 x 640 x 450 mm
<b>Weight</b>	Approx. 12.0 kg
<b>Power Supply</b>	AC 115 to 240 V
<b>Frequency</b>	50/60 Hz
<b>Power Consumption</b>	Less than 220 VA

### SPECIFICATIONS

<b>Stimulator Bowl Type</b>	Part hemispherical, radius 300 mm.
<b>Visual Field Extent</b>	Integrated diffusive screen 100°

#### Standard Test Fields

	green / white	blue / yellow
<b>Full 50°</b>	165 points	164 points
<b>Glaucoma 22°/50°</b>	101 points	100 points
<b>Central 30°</b>	117 points	116 points
<b>Central 22°</b>	93 points	92 points
<b>Wide 22°/30°</b>	128 points	127 points
<b>Peripheral 30° to 50°</b>	72 points	72 points
<b>Macula 10°</b>	45 points	44 points
<b>Driving 50°/80°</b>	189 points	188 points

<b>Stimulus Source</b>	Rear projection LED
<b>Stimulus Colour</b>	Green 570 nm, Blue 435 nm
<b>Stimulus Size</b>	Goldmann size III (0.43°) for green, Goldman size V for blue
<b>Stimulus Intensity</b>	0.03 asb to 1000 asb in 15 steps with 3 dB or 45 steps 1dB for green, 0.01 asb to 65 asb with 1 dB step for blue
<b>Exposure Time</b>	Adjustable: 0.1 to 9.9 s
<b>Response Time</b>	Adjustable: 0.1 to 9.9 s
<b>Inter Test Delay</b>	Adjustable: 0.1 to 9.9 s
<b>Background Illumination</b>	10 asb (3.2 cd/qm) for green on white test 100 cd/qm for blue on yellow test, automatic level control
<b>Fixation Control Method</b>	Heijl-Krakau - blind spot method and CCD camera
<b>Fixation Monitor</b>	Yellow LED, 588 nm
<b>Test Lens Diameter</b>	38 mm
<b>Chin Rest</b>	Electric
<b>PC</b>	Embedded, Pentium Celeron Mobile processor, 512 MB RAM, 40 GB HDD

#### Dimensions & Electric Requirements

<b>Dimensions WDH</b>	740 x 640 x 450 mm
<b>Weight</b>	Approx. 20.0 kg
<b>Power Supply</b>	AC 100 to 240 V
<b>Frequency</b>	50/60 Hz
<b>Power Consumption</b>	Less than 140 VA



**HS**  
**Optikmaschinen**

Handwerkerstraße 14  
48720 Rosendahl-Holtwick  
Tel: 02566/4720  
Fax: 02566/1620  
Email: [hsoptikmaschinen@hotmail.com](mailto:hsoptikmaschinen@hotmail.com)  
[www.hs-optikmaschinen.de](http://www.hs-optikmaschinen.de)